ID-Karte “DiaClon ABO/Rh for newborns”

A, B, AB, DVI+, ctl/ DAT

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ID-Card “DiaClon ABO/Rh for newborns”

A, B, AB, DVI+, ctl/ DAT

Determination of the ABO/Rh blood groups with direct antiglobulin test (DAT) for newborns

Product-Identification: 50961

| | |
|---------------------|---|
| Introduction | <p>Anti-A, anti-B and anti-AB test sera are necessary to detect the presence or absence of A/B antigens on human red cells. Since A and B antigens are not fully developed at birth, weaker reactions may occur with red cells of newborns than of adults and subgroups often cannot be identified.</p> <p>The serum of adults contains antibodies directed against the A and B antigens absent from their own red cells. Both antibodies appear after the first 4 to 6 months of life. As a result, reverse grouping is not usually undertaken on newborn blood samples. Confirmation of the newborn’s blood group is indicated when the A and B antigen expression is fully developed (2-4 years).</p> <p>The D antigen as well as weak D is fully developed at birth. The determination of the RhD status of the newborn’s blood group is important if the mother is RhD negative.</p> <p>A direct antiglobulin test (DAT) on newborn blood samples has become a standard procedure, since it of importance to know if the newborn’s red cells have been coated with maternal antibodies in-utero.</p> |
|---------------------|---|

| | |
|-----------------|--|
| Reagents | <div style="border: 1px solid black; padding: 2px; display: inline-block;">IVD</div> <p>ID-Card “DiaClon ABO/Rh for newborns” contains monoclonal anti-A [cell line LM 297/628 (LA-2)], anti-B [cell line LM 306/686 (LB-2)], anti-AB [cell lines ES131 (ES-15) + Birma-1 + ES-4] and anti-D (cell lines ESD-1M + 175-2), within the gel matrix. The microtube (ctl) is the negative control. The anti-human globulin serum is a blend of rabbit anti-IgG and monoclonal anti-C3d (cell line C139-9).</p> <p>Preservative: < 0,1% NaN₃.</p> <p><i>Caution: All reagents should be treated as potentially infectious.</i></p> <div style="display: flex; align-items: center;">  <div style="margin-left: 10px;"> <p><i>Do not store near any heat, air conditioning sources or ventilation outlets.</i></p> <p><i>Stability: see expiry date on label.</i></p> </div> </div> |
|-----------------|--|

| | |
|-------------------------------------|--|
| Additional reagents required | ID-Diluent 2: modified LISS for red cell suspensions. <i>(see related package insert)</i> |
|-------------------------------------|--|

| | |
|-----------------------------------|--|
| Further materials required | <ul style="list-style-type: none"> • ID-Dispenser • ID-Pipetor • ID-Tips (pipetor tips) • Tubes for suspensions • ID-Working table • ID-Centrifuge 6, 12 or 24 |
|-----------------------------------|--|

| | |
|------------------------|--|
| Sample material | <p>Draw blood samples using acceptable phlebotomy techniques. Cord or heel prick samples may be used. It is normally not necessary to wash the cells before use (see “Remarks”). Where cord samples are used, care must be taken to avoid contamination with Wharton’s jelly.</p> <p>Preferably, blood samples should be drawn into Citrate, EDTA or CPD-A. For reliable results, use of freshly collected blood is indicated.</p> |
|------------------------|--|

| | |
|------------------------------------|---|
| Preparation of blood sample | <p>Prepare a 0,8% red cell suspension in ID-Diluent 2 as follows:</p> <p><i>Allow the diluent to reach room temperature before use.</i></p> <ol style="list-style-type: none"> 1. Dispense 1,0 mL of ID-Diluent 2 into a clean tube. 2. Add 10 µL of packed red cells or 20 µL of whole blood, mix gently. <p><i>The cell suspension may be used immediately.</i></p> |
|------------------------------------|---|

| | |
|-----------------|---|
| Controls | Known positive and negative samples should be included in accordance with the relevant guidelines of quality assurance. |
|-----------------|---|

| | |
|-----------------------|---|
| Test procedure | <p>Do not use ID-Cards which show signs of drying, have bubbles, damaged seals, drops of gel or supernatant in the upper part of the microtubes or on the underside of the aluminium foil.</p> <ol style="list-style-type: none"> 1. Identify the ID-Card with the unique patient or donor number/details as appropriate. 2. Remove the aluminium foil from as many microtubes as required by holding the ID card in the upright position. 3. Pipette 50 µL of the red cell suspension to all 6 microtubes of the ID-Card. 4. Centrifuge the ID-Cards for 10 minutes in the ID-Centrifuge. 5. Read and record the results. |
|-----------------------|---|

Carte-ID “DiaClon ABO/Rh for newborns”

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Interpretation of results

A) Principle [2]

Positive: Agglutinated cells forming a red line on the surface of the gel or agglutinates dispersed in the gel.
 Negative: Compact button of cells on the bottom of the microtube.

B) Reactions for blood groups ABO

| Anti-A | Anti-B | Anti-AB | Blood group |
|-----------|-----------|-----------|-------------|
| + to ++++ | negative | + to ++++ | A |
| negative | + to ++++ | + to ++++ | B |
| + to ++++ | + to ++++ | + to ++++ | AB |
| negative | negative | negative | O |

In the presence of weak or very weakly expressed antigens the reaction can be negative. The anti-B of monoclonal origin does not react with the acquired B antigen.

Important: The microtube ctl must show a negative reaction. If the ctl is positive, the ABO determination is not valid. Repeat the test as described under "Remarks 1."

C) Reactions for RhD

| +++ to ++++ | ± to ++* | negative |
|--------------|-------------------|--------------|
| RhD positive | RhD weak positive | RhD negative |

* ±, trace or weak reactions should be subject to further investigations to distinguish between weak and partial D types as appropriate for the category of sample being tested.

The anti-D test sera were selected so as **to react** with DVI variants. Some DVI variants can give very weak reaction.

Important: The microtube ctl must show a negative reaction. If the ctl is positive, the RhD determination is not valid. Repeat the test as described under "Remarks 1."

D) Reactions for the direct antiglobulin test (DAT)

A negative reaction indicates the absence of detectable antibodies on the newborn's red cells.

A positive reaction (± to ++++) indicates that the newborn's red cells are sensitized (red cells coated with antibodies).

Remarks

- The negative control (ctl) must always show a negative reaction. If the negative ctl is positive, proceed as follows:
 - Wash the red cells first with isotonic saline solution (or ID-Diluent 2) before preparing the 0,8% cell suspension.
 - Then proceed as under "Preparation of blood sample" and "Test procedure".
- If the negative control remains positive, safe interpretation of the ABO and RhD groups cannot be made and further investigation is required.

Limitations

- ID-Cards which show air bubbles in the gel or drops in the upper part of the microtubes and/or the seal, must be centrifuged before use.
- Bacterial or other contamination of materials used can cause false positive or false negative results.
- Fibrin residues in the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
- Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
- Use of suspension solutions other than ID-Diluent 2 may modify the reactions.
- Too heavy or too weak red cell suspensions can cause aberrant results.

Bibliography

- Mollison, P.L., Engelfriet, C.P. and Contreras, M.: Blood Transfusion in Clinical Medicine, 10th ed. 1997; Blackwell Scientific Publications, Oxford.
- American Association of Blood Banks: Technical Manual; 11th ed. 1993.
- Lapierre, Y., Rigal, D., Adam, J. et al.: The gel test; A new way to detect red cell antigen-antibody reactions. Transfusion 1990; 30: 109–113.

Products

| | | |
|---------------------------------------|----------------|------------|
| ID-Card "DiaClon ABO/Rh for newborns" | 4 x 12 | REF 001047 |
| | 24 x 12 | REF 001048 |
| | 60 x 12 | REF 001049 |
| | 112 x 12 | REF 001050 |

These products are guaranteed to perform as described on the label and in the instruction sheet. The manufacturer declines all responsibility arising out of the use or sale of these products in any way or for any purpose other than those described therein.

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